

South Sound Radiology – Referring Physician Portal

Go to: <https://mobile.ssrportal.com> Enter your user name & password



ENTER YOUR PATIENT'S NAME *(it will display all patients with that name; last name, first name)*

- To select your patient, **Single** click anywhere on the study (or check the box to compare multiple studies / load data)
- The report and all available images for the study selected will be displayed in the right column
- **Double** click on the item you would like displayed

IMAGES

~ All images pertaining to the study you selected will be displayed on the left side; double click or drag-n-drop **the image or series you wish to view**

- If you would like to see all of the patient's studies, click on the icon in the bottom right corner of the column (hover over the icon – it's labeled **Related Studies**)

~ At the bottom of this column are several icons that will allow you to link and display images in various ways

~ Right click the image and a menu will pop up with other tools

~ To close your image and go back to search another patient, **DON'T USE YOUR BACK BUTTON**, click on the Study Browser button in the upper left corner

REPORTS

~ To print your report, click on the bottom right print icon

~ To close your report, click on the **X** in the top right corner

To open a User Guide, right click on any exam and select "**OPTIONS**", then click on the "**USER GUIDE**" button